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LAVANID® - Wundgel / LAVANID® - Wundgel V+

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is why treatment with preparations that are very well tolerated by the tissues is indicated.

LAVANID® Wound Gel and LAVANID® Wound Gel V+ are sterile, ready-to-use wound gels for cleaning and moistening wounds. LAVANID® Wound Gel V+ is more viscous than LAVANID® Wound Gel.


Polyhexanide 0,04% is used as a preservative. Polyhexanide is evaluated as the agent of choice for chronic wounds and burns. In addition to cleaning and moistening wounds, this wound gel is also appropriate for preservative dressings, e.g. for wounds colonised with multiresistant bacteria. Preservative dressings reduce the risk of carrying the bacteria from the wound to the patient's environment and they protect the wound from bacteria getting in from the outside.

SIDE EFFECTS

The ingredient polyhexanide can cause allergic reactions. There are individual case reports of acute immune system reactions (anaphylactic reactions). Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

INFORMATION ABOUT SHELF-LIFE

The expiry date can be found on the container. Do not use after the expiry date!

Shelf life after opening:

1 M 1 month

The wound gel is sterile as long as the packaging is unopened and undamaged. Although it is no longer sterile after the container has been opened, because it contains a preservative it can be used up to one month after opening under certain conditions. These conditions are that the contents are protected from contamination. The container must not come in contact with the wound. The opening should be capped immediately after each use. Gel residue around the opening should be wiped off with a sterile compress.

TISSUE TOLERANCE

Experimental and clinical studies have shown that the combination of Ringer's solution with 0.04% polyhexanide is very well tolerated by tissue.

Clinical trials with this wound gel have shown that it does not negatively affect the formation of healthy new tissue. There are no local signs of irritant conditions.

Use of this product is pain-free for patients.

Because it is so well tolerated by tissues the wound gel is suitable for repeated and long-term wound treatment.

DISPOSAL INFORMATION

Dispose of unused products or waste material in accordance with the specific national regulations.

INSTRUCTIONS FOR USE

- Prepare and clean the wound as usual.
- Apply the wound gel to the wound or a wound dressing (e.g. pads or compresses). Dress wound as usual. The wound gel should remain on the wound until the next dressing change. For optimal wound moistening we recommend increasing the amount of wound gel and lengthening the intervals between dressing changes. If the dressing sticks to the wound while being changed, moisten the dressing with a rinsing solution (e.g. LAVANID® Wound Rinsing Solution) before removing it.
- If the wound gel is going to be used more than once, please read the shelf-life information.

CONTRAINDICATIONS

The wound gel must not be used:

- in cases of known hypersensitivity to any of the ingredients,
- in the area of cartilage and joints.

There is insufficient data on its use during pregnancy and breast-feeding. Therefore, this wound gel should not be used unless urgently indicated.

WARNINGS AND PRECAUTIONS FOR USE

For external use only.
Do not mix with additives.
Sterile if package is undamaged.
Do not sterilize! This may critically alter its usability.
Do not use after the expiry date.

DESCRIPTION AND PRODUCT PROFILE

The surface of chronic wounds is generally made up of devitalized and necrotic tissue, and is also colonised with bacteria. For wound treatment to be successful the wound surface needs to be removed, e.g. through surgical débridement. Wound irrigation and moistening are a necessary part of this process. Chronic wounds exhibit dysfunctional wound healing processes. This

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