


Kurzzugbinde
Bande à allongement court
Short Stretch Bandage
Venda de tracción corta
Ligadura de baixa elasticidade
Benda ad allungamento limitato
Korte-rekzwachtel
Kortstræksbind
Lågælastisk binda
Obinadlo s krátkým tahom
Ovinadlo s krátkým tahom
Rövid megnyúlású pólya
Bandaž o krótkim naciągu
Kısa cekisli bandaj

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Durelast® / Durelast® Combi: 66% cotton, 34% polyamide

Intended use

For compression, retention of dressings, support and relief and for immobilisation of limbs.

Indications

phlebological and lymphological indications: varicosis, venous insufficiency, thrombophlebitis, phlebothrombosis, post-thrombotic syndrome, stage I to III b chronic venous insufficiency according to Widmer, venous and mixed ulcers, oedema (phleboedema, lymphoedema, lipoedema, post-traumatic, post-operative, during pregnancy), after sclerotherapy and venous surgery. contusions, sprains, tendon injuries and joint dislocations

Contraindications

Advanced peripheral arterial occlusive diseases, congestive heart failure, septic phlebitis, phlegmasia cerulea dolens, sensation disorders of the skin, known allergies and/ or hypersensitivity to any of the product components.

Side effects

Skin irritations and/or allergies can occur in very rare cases.

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Warnings and precautions

Incorrect bandaging can result in skin necrosis and pressure damage to nerves. Slight bluish tinge of the toes is normal, this should disappear as soon as the patient starts to walk. If the symptoms do not disappear (pain and bluish tinge), the bandage must be removed and applied again with less compression. The patient must inform their doctor if sharply increasing pain occurs while wearing the bandage. For use with ABPIs over 0.8. When ABPI is between 0.5 and 0.8, apply under medical supervision only. Use adhesive strips to fix the bandage in place. Do not fix in place with bandage clips because these can cause injury.

General instructions

- Before applying the bandage, clean/debride any wound that is present and cover with a wound dressing appropriate for the wound status.
- The bandaging technique is chosen in consultation with the treating physician.
- It is recommended to apply the bandage in the morning before the patient gets up or after the patient has rested the leg/arm in a raised position for a prolonged period and to then encourage

the patient to walk/move their arm.

- The width of the bandage should correspond to the diameter of the limb to be bandaged.
- To achieve maximum joint movement, it is recommended to apply the bandage with the foot at an angle of 90° to the leg.
- The bandage must be applied so that the pressure from the bandage decreases from the foot to the knee and, where applicable, to the gluteal fold.
- The bandage is unrolled with the rolled-up portion on top. Hold the top of the bandage close to the leg/arm.

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bandage may increase with washing. Bear this in mind when applying the bandage. Durelast®: Can be washed up to 10 times at 60°C (standard commercially available mild detergent/delicates cycle). Not suitable for tumble drying. Do not iron.

Disposal

In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance – AVV). Packaging labelled as recyclable should be taken to the applicable national recycling systems.

General instructions

If a serious incident occurs, contact the manufacturer and the respective health authorities.

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