Idealbinde

Idealbinde steril

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en Instructions for Use

Product description and performance characteristics:

- Idealbinde: Stretch approx. 90 %
- Lenkideal®: Stretch approx. 85 %

Product composition:

- Idealbinde: 100% cotton
- Lenkideal®: 64% cotton, 34% polyamide, 2% elastane

Intended purpose:

For fixation, support and relief, for immobilisation of limbs and for compression.

Indications:

phlebological and lymphological indications: varicosis, deep venous insufficiency, thrombophlebitis, phlebothrombosis, post-thrombotic syndrome, Widmer's stage I to III b chronic venous insufficiency, venous and mixed ulcers, oedema (phleboedema, lymphoedema, lipoedema, post-traumatic, post-operative, during pregnancy), after sclerotherapy and venous surgery. contusions, sprains, tendon injuries and joint dislocations

Contraindications:

Advanced peripheral arterial occlusive diseases, congestive heart failure, septic phlebitis, phlegmasia cerulea dolens, sensation disorders of the skin, known allergies and/or hypersensitivity to any of the product components.

Side effects:

In very rare cases skin irritations and/or allergies may occur.

Warnings and precautions:

- Incorrect bandaging can result in skin necrosis and pressure damage to nerves.
- Slight bluish tinge to the toes is normal, this should disappear as soon as the patient starts to walk.
- If the symptoms do not disappear (pain and bluish tinge), the bandage must be removed and applied again with less compression.
- The patient must inform their doctor if sharply increasing pain occurs while wearing the bandage.
- For use with ABPIs over 0.8. When ABPI is between 0.5 and 0.8, apply under medical supervision only.
- Use adhesive strips to fix the bandage in place. Do not fix in place with bandage clips because these can cause injury.
- If appropriate, an underlying bandage is advisable – depending on the condition of the skin.
- Known allergy and / or hypersensitivity to any of the product components.

General instructions:

- Before applying the bandage, clean/debride any wound that is present and cover with a wound dressing appropriate for the wound status.
- It is recommended to apply the bandage in the morning before the patient gets up or after the patient has rested the leg in a raised position for a prolonged period and to then encourage the patient to walk.
- The width of the bandage should correspond to the diameter of the limb to be banded.
- To achieve maximum joint movement, it is recommended to apply the bandage with the foot at an angle of 90° to the leg.
- The bandage must be applied so that the pressure from the bandage decreases from the foot to the knee and, where applicable, to the gluteal fold.
- The bandage is unrolled with the rolled-up portion on top. Keep the roll close to the leg/arm.

Usage information:

- Idealbinde and Lenkideal® can be applied using different compression techniques (illustrations A1-A11 or B1-B10).
- The bandaging technique is chosen in consultation with the treating physician.

Care recommendation:

Washing may cause the ends of the bandage to fray. It will retain its short-stretch characteristics and function. The maximum stretch of the bandage may increase with washing. Bear this in mind when applying the bandage. Idealbinde: Can be washed up to 5 times at 95oC (standard commercially available mild detergent/delicates cycle).

Disposal:

In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance - AVV). Packaging labelled as recyclable should be taken into the respective national recycling systems.

General instructions:

In case of a serious incident, contact the manufacturer and the responsible health authorities.

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