
Komprex®

Schaumgummi zur Polsterung bei Kompressionsverbänden
Caoutchouc-mousse de rembourrage pour bandages compressifs
Foam Rubber for padding under compression bandages
Gomaespuma para almohadillado en vendas de compresión
Espuma de borracha para acolchoamento de ligaduras de compressão
Gommapiuma per imbottitura di bendaggi compressivi
Schuimrubber voor polstering bij compressieverbanden
Skumgummi til polstring ved kompressionsforbindinger
Skumgummi som polstring till kompressionsförband
Pěnová guma k polstrování při kompresních obvazech
Penová guma na vypolstrovanie pri kompresných obvázoch
Habgumi kompressziós kötések kipárnázásához
Guma piankowa do wyściełania opatrunków uciskowych

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en Instructions for use

Product description and performance characteristics:

Komprex® Foam Rubber Pads: Ready-to-use
Komprex® Foam Rubber Bandages: Ready-to-use and can be cut to size
Komprex® Foam Sheets: can be cut to size

Product composition:

Komprex® Foam Rubber Pads: styrene-butadiene rubber
Komprex® Foam Rubber Bandages: 80% styrene-butadiene rubber
Komprex® Foam Sheets: 80% styrene-butadiene rubber

Intended purpose:

For padding under compression bandages

Indications:

Phlebological and lymphological indications: Varicosis, deep venous insufficiency, thrombophlebitis, phlebotrombosis, post-thrombotic syndrome, Widmer’s stage I to III b chronic venous insufficiency, venous and mixed ulcers, oedema (phleboedema, lymphoedema, lipoedema, post-traumatic, post-operative, during pregnancy), after sclerotherapy and venous surgery. Contusions, sprains, tendon injuries and joint dislocations.

Contraindications:

Advanced peripheral arterial occlusive diseases, congestive heart failure, septic phlebitis, phlegmasia cerulea dolens, sensation disorders of the skin, known allergies and/or hypersensitivity to any of the product components.

Side effects:

In (very) rare cases, skin irritations and/or allergies may occur.

Warnings and precautions:

- Incorrect bandaging can result in skin necrosis and pressure damage to nerves.
- Before applying the bandage, clean/debride any wound that is present and cover with a wound dressing appropriate for the wound status.
- Not intended to be sterilised.

Usage information:

- A. Komprex® Foam Rubber Pad:
1. a) Wrap Komprex® Pad with gauze
 - b) Wound pad
 - c) Ulcer
2. Apply the rounded side in direction to the body without tension
3. Apply the compression bandage using an appropriate wrapping technique
- B. Komprex® Foam Rubber Bandage:
1. Wrap the limb with a bandage liner

- Wrap Komprex® Bandage without pulling tight or overlapping or cut to size
2. Apply the compression bandage using an appropriate wrapping technique
- C. Komprex® Foam Rubber Sheet:
1. Cut out the desired shape
 2. Angle the edges
 3. For the remaining application steps, see Komprex® Foam Rubber Pad

Disposal:

In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance - AVV). Packaging labelled as recyclable should be taken into the respective national recycling systems.

General instructions:

In case of a serious incident, contact the manufacturer and the responsible health authorities.

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