


Raucocel®

Nasentamponade
Méchage Nasal
Nasal Packing
Taponamiento nasal
Tamponamento nasal
Tampone nasale
Neustampon
Nasetampon
Nästamponad
Nosní tamponáda
Nosová tamponáda
Orrtampon
Tamponada za nos
Pivukó tqurúv
Ninaverejooksu tampoon

Huizhou Medical Devices Co., Ltd.
North Shangxia Rd. Dongjiang
Hi-tech Industry Park, 516005,
Huizhou, P.R.China

Shanghai International Holding Corp.
GmbH (Europe), Eiffelstrasse 80,
20537 Hamburg, Germany

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(or) near fainting when standing up), and/or a rash looking like a sunburn. Nasal packing exhibit no antimicrobial properties. Rarely, displacement and ingestion of a nasal packing may occur. The expansive sponge is non-toxic and is considered safe to pass once it is completely in the alimentary tract. The patient should be monitored in the normal manner for the ingestion of any non-toxic material. If aspiration should occur, prompt medical intervention is imperative with close patient monitoring and airway support pending the removal of the material.	
Warnings and precautions Packaging must be checked before usage. Do not use if the packaging is damaged. Keep it dry once opened. The product should not be moistened with water or water-bearing liquid before use. When the product is used for ENT department, the string must be fixed to the wing of the nose or the cheek to prevent being aspirated into the respiratory airway. Do not keep on the wound for longer than 72 hours. This product is not reusable.	
Application 1. Select a proper type of Raucocel® Nasal Packing and unpack the product.	
Indications Raucocel® Nasal Packing for post-operative hemostasis of nasal cavity and paranasal sinus.	
Contraindications Patients with allergies to any of the components should refrain from using this product.	
Adverse effects In rare instances the physiochemical conditions associated with nasal surgery, both with and without nasal packings, may present a risk of toxic shock syndrome (TSS). Warning signs of TSS include sudden fever, vomiting, diarrhea, dizziness, fainting	
Product composition Polyvinyl acetal, string, airway tube: silicone rubber.	

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