Rosidal® TCS

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en Instructions for Use

Product description and performance characteristics

Two-component compression system

- **Rosidal® SC:** padding compression bandage (stretch approx. 50%; white), which is applied directly to the skin.
- **Rosidal® CC:** natural latex-free, cohesive compression bandage (stretch approx. 75%; cream-coloured) with short stretch for therapeutically effective compression.

Product composition

- **Rosidal® SC:** Polyurethane foam with textile fabric made of 60% polyamide, 40% cotton
- **Rosidal® CC:** 51% cotton, 39% polyamide, 10% polyurethane, with a synthetic rubber coating

Intended purpose

For compression, retention of dressings, support and relief and for immobilisation of limbs.

Indications

Venous insufficiency according to Widmer, venous and mixed ulcers (ankle-brachial index (ABPI) above 0.8 or with an ABPI of 0.5-0.8 under medical control by a physician), chronic oedema.

Contraindications

Advanced peripheral arterial occlusive

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diseases, decompensated cardiac insufficiency, septic phlebitis, phlegmasia coerulea dolens, perception disorders of the skin, intolerance to the materials of the bandages.

Side effects

In very rare cases skin irritations and/or allergies may occur.

Warnings and precautions

- Cannot be used on ankle circumferences smaller than 18 cm, unless the padding material increases the circumference to over 18 cm.
- Incorrect bandaging can result in skin necrosis and pressure damage to nerves.
- Slight bluish tinge to the toes is normal, this must disappear as soon as the patient starts to walk.
- If the symptoms do not disappear (pain and bluish tinge), the bandage must be removed and applied again with less compression.
- The patient must inform their doctor if sharply increasing pain occurs while wearing the bandage.
- In the case of diabetes with advanced microangiopathy – only use under medical supervision.
- Do not fix in place with bandage clips because these can cause injury.

General instructions

It is recommended to apply the bandage in the morning before the patient gets up or after the patient has rested the leg in a raised position for a prolonged period and to then encourage the patient to walk.

- The bandage must be applied so that the pressure from the bandage decreases from the foot to the knee.
- The bandage can be worn during the day and at night.
- The compression bandage can stay on the leg for up to 7 days, depending on the wound status and wound dressing changes.
- The bandage can be removed by unwrapping it or simply cutting it off.

Usage information

Préparation

- The usual diagnostic examinations to determine the vascular status must be carried out before applying the bandage.
- Before applying the bandage, clean/ debride any wound that is present and cover with a wound dressing appropriate for the wound status.

Application instructions

- Training by healthcare professionals on the correct application is recommended.
- The bandage is unrolled with the rolled-up portion (roll) on top. Keep the roll close to the leg.
- Apply both bandages at full stretch (max. stretch = Safe-Loc®).
- To achieve maximum joint movement, it is recommended to apply the bandage with the foot at a right angle (90°) to the leg (Fig. 1).
- Excess bandage material can be cut off.

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- Use adhesive strips to fix the bandage in place. Do not fix in place with bandage clips because these can cause injury.

- a. Measure the ankle circumference. If necessary, increase the circumference of the leg (if <18 cm) with the padding compression bandage (white).
- b. Wrap the padding compression bandage at full stretch twice around tarsus (Fig. 2), including the heel (Fig. 3). Continue, overlapping by approx. 50% until about 2 cm below the head of the fibula. Affix the end (Fig. 4).
- c. Wrap the cohesive compression bandage (cream-coloured) exactly like the first layer, first twice around tarsus at full stretch (Fig. 5).
- d. Roll the bandage using a figure-of-eight technique around the back of the leg, including the heel (Fig. 6) until about 3 cm below the head of the fibula.
- e. Finished bandage (Fig. 7).
- f. To prevent the edges of the bandage from rolling up, they can be fixed in place with long adhesive bandage strips.

Disposal

In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance – AVV). Packaging labelled as recyclable should be taken into the respective national recycling systems.

General instructions

In case of a serious incident, contact the manufacturer and the responsible health authorities.

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