

nitrile examination gloves ultra comfort

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ENGLISH

Product Description

Klinion Powder-Free Nitrile Examination Gloves are disposable gloves made of nitrile butadiene rubber that is intended to wear on the hand for medical purposes to provide a barrier against potentially infectious materials and other contaminants. These gloves are also intended to protect against risks associated with contact against certain chemicals, chemotherapy drugs, and microorganisms (bacteria, fungi and virus), where hand protection is of priority. Not made with natural rubber latex.

Directions for use

- Gloves are packed in dispenser which is suitable for transport. Keep the gloves in the box when not in use and keep away from fluorescent- and sunlight.
- Before usage, inspect the gloves for any defects or imperfections.
- Usage: For single use only. If re-used: Extremely high risk of cross-contamination, deterioration of barrier protection, deterioration of glove's properties and loss of lot traceability.
- Sizing – Select the right size glove for your hand.
- Donning – Hold glove by the bead with one hand. Align the glove thumb with your other hand thumb and slide your hand into the glove, one finger into each glove finger. Pull by the glove palm to a get a good fit. Don the other glove by the same procedure.
- Inspection – Punctures or tears may occur after donning. Inspect each glove after donning, and immediately discontinue use if found damaged.
- Doffing – Hold glove bead and pull toward the finger until the glove comes off.
- Disposal – Properly disposal of all used gloves. Follow the local regulations and your Institution's policies for disposal.

Performance and Limitation of Use

This product had been tested in accordance with EN ISO 374-1:2016 + A1:2018, EN ISO 374-2:2019, EN 374-4:2019, EN 16523-1:2015+A1:2018. This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.

This product had been tested in accordance with EN ISO 374-5:2016, and ISO 16604:2004 to meet the requirements for resistance to penetration by blood-borne pathogens, test method using Phi-X174 bacteriophage. This product provides protection against bacteria, fungi and virus. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

The gloves were found to meet with the REACH annex XVII requirements for Polycyclic Aromatic Hydrocarbons (PAHs).

Components used in glove manufacturing may cause allergic reactions to some users. If allergic reactions occur, seek for medical advice immediately. Where relevant, a list of the substances contained in the glove which are known to cause allergies, per listed in Annex G of EN ISO 21420:2020, shall be supplied on request.

In case a serious incident has occurred in relation to the device, inform Medeco BV as listed on the packaging, your local distributor and your local healthcare authority.

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