

# UROMED SUPRAPUBISCHER BALLON-KATHETER

## SUPRAPUBIC BALLOON CATHETER

### REF 2956/2966

Kurze, zentral offene Spitze  
2 Augen


*Short, central open Tip  
2 eyes*

### Ballon-Katheter zur suprapubischen Blasendrainage

- kurze, atraumatische Spitze mit abgerundeter Kante
- flache Ballonform
- großlumig
- Soft-Ballon
- Silikon
- weich

### Ballon Catheter for suprapubic drainage of the urinary bladder

- *short, atraumatic tip with rounded edge*
- *flat shape of balloon*
- *wide lumen*
- *soft balloon*
- *silicone*
- *soft*



**REF 2956**  
Länge 43 cm  
*length 43 cm*

**REF 2966**  
Länge 33 cm  
*length 33 cm*

REF <i>REF</i>	Größen <i>Sizes</i>	Ballon <i>Balloon</i>	PZN <i>PZN</i>
2956-10	Ch. 10	2 ml	00051724
2956-12	Ch. 12	5 ml	00051730
2956-14	Ch. 14	5 ml	00051747
2956-16	Ch. 16	5 ml	00051753
2956-18	Ch. 18	5 ml	00051776
2956-20	Ch. 20	5 ml	00051807
2956-22	Ch. 22	5 ml	02784388
2956-24	Ch. 24	5 ml	02784394

REF <i>REF</i>	Größen <i>Sizes</i>	Ballon <i>Balloon</i>	PZN <i>PZN</i>
2966-10	Ch. 10	2 ml	07654443
2966-12	Ch. 12	5 ml	04741957
2966-14	Ch. 14	5 ml	03281465
2966-16	Ch. 16	5 ml	04741963
2966-18	Ch. 18	5 ml	04741986
2966-20	Ch. 20	5 ml	04741992
2966-22	Ch. 22	5 ml	02784402
2966-24	Ch. 24	5 ml	02784431

**Verpackung:** doppelt, steril  
**Original-Karton:** 10 Stück

**Packing:** double, sterile  
**Original carton:** 10 pieces

**UROMED**  
PRODUKTE FÜR DIE UROLOGIE