



**NOBAMED Paul Danz AG**

## NOBAHEBAN<sup>®</sup>

**REF 455405**

### **Product Description, Intended use, Application**

The blue cohesive compression bandage (size 5 cm x 4.5 m) has a stretchability of approx. 100 %. It is used for bandages for short-term compression. It is applied over a gauze ball or swab. It is also used as a supporting and relieving bandage

### **Composition**

Polyamide, polyester, latexfree elastomere

### **Contraindications**

The product should not be used in the case of a known allergy against the material.

### **Note**

During application, cohesive bandages should not be stretched too much, as otherwise the arterial blood supply might be disturbed.

### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### **Normative and Regulative Requirements, Common Standards**

Medical Device according to *MDD 93/42/EEC, MDR (EU) 2017/745*.

The product does not contain dangerous toxic substances according to REACH.

### **Packaging**

Primary: foil

Secondary: folding box  
made of cellulose

Tertiary: carton made of  
cellulose

### **Symbols used in labelling**

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

### **Storage and Transport**

Dry and dustfree

### **Single use device**

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis.